


Deluxe Viral Disinfection Kit



Inactivates Viral Pathogens

This kit is designed to allow anyone to quickly make disinfectant solutions which are effective against contagious viral and bacterial pathogens including Human Corona Virus, Influenza, Adenovirus and many others.

64 Millennium Q Disinfectant is a hospital grade disinfectant-cleaner that when used properly, will inactivate Human Corona Viruses, including the novel human corona virus identified as SARS-CoV-2, the cause of COVID-19, when used in accordance with label instructions for disinfecting against Human Corona Virus.

Check List of High Touch Critical Surfaces:

- Bedrails
- Tables
- Light switches
- Door knobs
- Hand rails
- Window sills
- Desk tops
- Chairs
- Telephones
- Door frames



Makes 4 full quarts of Hospital Grade Disinfectant.

The Kit includes:

- 1 Viral disinfectant concentrate
- 1 Mixing Bottle
- 4 Empty quart bottles
- 4 Trigger sprayers
- 4 Microfiber Towels
- Detailed Instructions



Fill mixing bottle with water up to within 2.5 inches from top. Empty the 2 oz bottle of Millennium Q into the mixing bottle, place cap on bottle, gently mix.



Fill quart bottles with ready to use solution you just made. Attach trigger sprayer.



Spray all hard surfaces, allow to remain wet for ten minutes.



600 Cardigan Road
Shoreview, Minnesota
Phone 651-481-1900
Web www.multi-clean.com

Order Part number 903020

21032020